

# HEMOCHRON Signature Elite

Point-of-Care Automated  
Coagulation Analyzer



## FEATURES:

- ▶ Measured items: ACT, aPTT, PT/INR
- ▶ Sample volume: <50 µl whole blood
- ▶ Results are displayed on the screen and/ or printed
- ▶ ACT+ results are NOT influenced by hemodilution or hypothermia
- ▶ Portable

**SIGNATURE ELITE**  
**Accriva**  
diagnostics

## CONSUMABLES:

Product name	Notes
Cuvettes for the determination of ACT – ACT+	Sensitive to heparin levels of 1 to 6 U/ml. Ideal for cardiac surgery centers
Cuvettes for the determination of ACT – ACT-LR	Sensitive to heparin levels of 0 to 2.5 U/ml Ideal for cardiac catheterization centers, hemodialysis, hemofiltration and ECMO
Cuvettes for the determination of PT/INR	Designed for monitoring the oral anticoagulation therapy. Ideal for the Emergency Room, AIC, the cardiology and neurology departments
Cuvettes for the determination of PT/INR in citrate blood samples	
Cuvettes for the determination of aPTT	Sensitive to heparin levels of 0 to 1.5 U/ml. Ideal for the Emergency Room, AIC, the cardiology and neurology departments
Cuvettes for the determination of aPTT in citrate blood samples	

**Worldwide there are more Hemochron instruments in hospitals than any other type of POC medical equipment.**

**SYNTERGY CONSULT**



**BUCUREȘTI:** Strada Radovanu nr 31-33, sector 2,  
tel: +40 21 252 23 58/ 59, Fax : +40 21 252 07 60

**CLUJ-NAPOCA:** Aleea Prof. Valeriu Bologa Nr. 3, Bl.2, Parter,  
Telefon: +40 264 594 517, Fax: +40 264 447 507

www.synttergy.ro  
office@synttergy.ro